

Electric wheelchair G02

Product specification	
Material	Steel
Driving range	15~20km
Max speed	6km/h
Overall size	1095*630*965mm
Fold size	770*370*760mm
Seat size	450*430mm
Loading capacity	100kg
Front tire	10 inch (10"×2.00) ---solid
Rear tire	18 inch (18"×2.125") ---air
Controller	Domestic
Motor	DC 250W *2Pcs spur gear motor
Parking brake	Manual brake
Battery choice	24V 12Ah lead acid
N.W./G.W.	44kg/50kg
Common description:	
<p>Larger front wheel, better obstacle climbing capacity</p> <ol style="list-style-type: none"> 1. Foldable and portable, easy fold and put into car boot 2. Rear wheel drive by two motors 3. Easy to operate, dismantle and assembling 4. Removable battery box for convenient recharging 5. Flip-up armrest, adjustable and detachable footrest 6. Controller can be exchanged from right side to left side for both hands use 	
Detailed description /advantages:	
Frame	Full set of robot-welded steel frame. Liquid painting, smooth and shining
Motor	Domestic spur gear drive motor without electromagnetic brake
Battery	Chinese first class brand battery ---Chilwee. Back-mounted design, easy to remove.
Brake	Manual brakes on two sides, easy operate
Wheel	Front wheel: 10" (10"×2.00) high quality PU wearable solid tire Rear wheel: 18" (18"×2.125") high quality wide air tire
Backrest	Foldable back, plastic coated trigger
Armrest	Flip-up armrest, PU arm pad, desk armrest design
Footrest	Adjustable and detachable footrest
Cushion	Removable and washable cushion
Safe belt	Enhance the using safety

